

	MedicVOD IHE Integration Statement		
	PSW-81-001		
			Rev.: 01

MedicVOD

IHE Integration Statement

IHE Integration Statement			
Vendor	Product Name	Version	Date
Orpheus Medical	MedicVOD [MVD-0001]	01	October 2015
This Product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer	None	
	Imaging Document Source	None	
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer	None	
Internet address for vendor's IHE Information: www.orpheus-medical.com			
Links to Standards Conformance Statements for the Implementation			
HL7	www.orpheus-medical.com/ihe		
DICOM	www.orpheus-medical.com/ihe		
Links to general Information on IHE			
In North America: www.ihe.net	In Europe: http://www.ihe-europe.net/	In Japan: www.ihe-i.org	